


Dersalazine sodium

Chemical Properties

CAS No. :	367249-56-7
Formula:	C35H32N6NaO4
Molecular Weight:	623.669
Appearance:	no data available
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year



Biological Description

Description	Dersalazine, a PAF antagonist and TNF α inhibitor, is used potentially for the treatment of ulcerative colitis. It showed anti-inflammatory effect in reactivated colitis in rats, as evidenced by a significant reduction in MPO activity.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6034 mL	8.0171 mL	16.0341 mL
5 mM	0.3207 mL	1.6034 mL	3.2068 mL
10 mM	0.1603 mL	0.8017 mL	1.6034 mL
50 mM	0.0321 mL	0.1603 mL	0.3207 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481