

BCAT-IN-2

Chemical Properties

CAS No. : 1800024-45-6


Formula: C₁₇H₁₄ClF₂N₅O

Molecular Weight: 377.776

Appearance: no data available

Storage: store at low temperature

Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	BCAT-IN-2 is an orally active and selective inhibitor of mitochondrial branched-chain aminotransferase (BCATm), inhibits BCATm and BCATc, and is used in obesity and dyslipidemia.
Targets(IC50)	Others
In vitro	BCAT-IN-2 inhibits BCATm in primary human adipocytes, with a pIC ₅₀ of 6.5[1].
In vivo	BCAT-IN-2 (5 mg/kg, po; 1 mg/kg, iv) exhibits favorable pharmacokinetic parameters in mice[1].

Solubility Information

Solubility	DMSO: 30 mg/mL (79.41 mM), Sonification is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.647 mL	13.2352 mL	26.4704 mL
5 mM	0.5294 mL	2.647 mL	5.2941 mL
10 mM	0.2647 mL	1.3235 mL	2.647 mL
50 mM	0.0529 mL	0.2647 mL	0.5294 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Bertrand SM, et al. The Discovery of in Vivo Active Mitochondrial Branched-Chain Aminotransferase (BCATm) Inhibitors by Hybridizing Fragment and HTS Hits. J Med Chem. 2015 Sep 24;58(18):7140-63.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481