

Sprengerinin A

Chemical Properties


CAS No. : 88866-99-3

Formula: C₃₈H₆₀O₁₂

Molecular Weight: 708.88

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Sprengerinin A is a natural product of Ophiopogon, Liliaceae. The catalog number is TN5296 and the CAS number is 88866-99-3. Sprengerinin A can be used as a reference standard.
-------------	----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.4107 mL	7.0534 mL	14.1068 mL
5 mM	0.2821 mL	1.4107 mL	2.8214 mL
10 mM	0.1411 mL	0.7053 mL	1.4107 mL
50 mM	0.0282 mL	0.1411 mL	0.2821 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:36 Washington Street,Wellesley Hills,MA 02481