

exo-BCN-PEG2-Boc-amine

Chemical Properties


CAS No. : 1807501-87-6

Formula: C₂₂H₃₆N₂O₆

Molecular Weight: 424.53

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Exo-BCN-PEG2-Boc-amine is a PEG derivative containing a BCN group and a Boc-protected amine. PEG Linkers and PEG derivatives may be useful in the development of antibody-drug conjugates.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3555 mL	11.7777 mL	23.5555 mL
5 mM	0.4711 mL	2.3555 mL	4.7111 mL
10 mM	0.2356 mL	1.1778 mL	2.3555 mL
50 mM	0.0471 mL	0.2356 mL	0.4711 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Sano K, Nakajima T, Miyazaki K, Ohuchi Y, Ikegami T, Choyke PL, Kobayashi H. Short PEG-linkers improve the performance of targeted, activatable monoclonal antibody-indocyanine green optical imaging probes. *Bioconjug Chem.* 2013 May 15;24(5):811-6. doi: 10.1021/bc400050k. Epub 2013 May 3. PubMed PMID: 23600922; PubMed Central PMCID: PMC3674550.

Harrison E, Coulter JA, Dixon D. Gold nanoparticle surface functionalization: mixed monolayer versus hetero bifunctional peg linker. *Nanomedicine (Lond).* 2016 Apr;11(7):851-65. Review. PubMed PMID: 27021417.

Augusto MT, Hollmann A, Porotto M, Moscona A, Santos NC. Antiviral Lipopeptide-Cell Membrane Interaction Is Influenced by PEG Linker Length. *Molecules.* 2017 Jul 15;22(7). pii: E1190. doi: 10.3390/molecules22071190. PubMed PMID: 28714870; PubMed Central PMCID: PMC5776016.

Tuma R, Russell M, Rosendahl M, Thomas GJ Jr. Solution conformation of the extracellular domain of the human tumor necrosis factor receptor probed by Raman and UV-resonance Raman spectroscopy: structural effects of an engineered PEG linker. *Biochemistry.* 1995 Nov 21;34(46):15150-6. PubMed PMID: 7578129.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481