

(+)- Praeruptorin C

Chemical Properties


CAS No. :

Formula: C₂₂H₂₂O₇

Molecular Weight: 398.4

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	(+)- Praeruptorin C extracted from <i>Peucedanum praeruptorum</i> Dunn roots.
Targets(IC ₅₀)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.510 mL	12.5502 mL	25.1004 mL
5 mM	0.502 mL	2.510 mL	5.0201 mL
10 mM	0.251 mL	1.255 mL	2.510 mL
50 mM	0.0502 mL	0.251 mL	0.502 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:36 Washington Street,Wellesley Hills,MA 02481