


NU-7031

## Chemical Properties

CAS No. : 79105-88-7  
Formula: C<sub>14</sub>H<sub>15</sub>NO<sub>4</sub>  
Molecular Weight: 261.27  
Appearance: no data available  
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



## Biological Description

Description	NU-7031 is a potent DNA-PK inhibitor with potential anticancer activity.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.8275 mL	19.1373 mL	38.2746 mL
5 mM	0.7655 mL	3.8275 mL	7.6549 mL
10 mM	0.3827 mL	1.9137 mL	3.8275 mL
50 mM	0.0765 mL	0.3827 mL	0.7655 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481